

LAB ARMOR® BEAD BLOCKS™



The next generation of dry baths. In full color.

We are replacing old things in the lab again. The new Bead Block is no ordinary dry bath heat block. It's not designed around a solid block, but around Lab Armor beads for better flexibility and performance. Available in both single and double block sizes.



Portable



Fits different sized samples



Choose your color

Single Bead Block



Technical Specifications

Length: 3.75 inches (9.5 cm)

Width: 3.0 inches (7.6 cm)

Height: 3.0 inches (7.6 cm)

Weight: 2 lbs (0.95 kg)

Tub (LWH): 3.25 x 2.5 x 2.5 inches
(8.4 x 6.4 x 6.4 cm)

Beads: 0.25 L

Double Bead Block



Technical Specifications

Length: 3.75 inches (9.5 cm)

Width: 6 inches (15.2 cm)

Height: 3 inches (7.6 cm)

Weight: 4 lbs (1.9 kg)

Tub (LWH): 3.25 x 5.5 x 2.5 inches
(8.4 x 14 x 6.4 cm)

Beads: 0.50 L

Simple & Versatile



Thermal Uniformity:

At 95°C +/- 1°C (with lid)

Temperature Range:

From -80°C to 180°C

Too many blocks?

With a Lab Armor Bead Block, one size fits all.

Bead Block™ by Lab Armor is an aluminum block filled with Lab Armor Beads™ designed for use with any standard commercially available dry bath. A Bead Block replaces common drilled-out aluminum heat blocks and provides for optimal temperature uniformity and accuracy while incubating common lab vessels, including microfuge tubes, test tubes, and even microplates and slides.

Lab Armor Single Bead Blocks

Blue: 52100-BLU
Gold: 52100-GLD
Black: 52100-BLK
Silver: 52100-SLV
Red: 52100-RED

Lab Armor Double Bead Blocks

Blue: 52200-BLU
Gold: 52200-GLD
Black: 52200-BLK
Silver: 52200-SLV
Red: 52200-RED



labarmor.com | sales@labarmor.com | +1-469-496-2350